


Flocare® Feeding Set

Flocare® feeding sets are designed to safely and securely connect feeds to feeding tubes

Features

- The Flocare® pack feeding set connects to Nutrison and Nutrini packs by a plus shaped spike end, which can only be connected to enteral feed packs and containers, increasing patient safety
- Flocare® universal feeding sets connect to bottles by a screw connector/adaptor
- All Flocare® feeding sets are compatible with the Flocare® Container
- Accepts a female luer syringe
- Roller clamp
- Latex and DEHP free materials
- Coloured lilac to make a clear distinction from parenteral feeding systems



Recommendations for Use

- Flocare® feeding sets are intended for single use only and should be changed every 24 hours

Ordering information

Product	Nappi Code	Units per carton
Flocare Infinity Universal Pack Set with adaptor	683223001	30

Date of publication: May 2016

A food for special medical purposes; to be used under strict medical supervision.

Nutricia Southern Africa (PTY) Ltd
199 Bryanston Drive, Bryanston, 2191

For more information call the
Nutricia Customer Care Line 087 350 2310

